



This versatile ChemInstruments software operates multiple test platforms with point and click commands and provides test data results in report ready format. With the power of your computer, EZ-Lab software brings test results directly into your harddrive providing powerful statistical analysis and report generation.

- One common screen contains all commands, including the calibration sequence.
- Test data is graphed immediately upon test completion.
- The graph cropping feature provides more detail



analysis and selected views.

- Files can be saved to two additional locations.
- Test data files can be transferred across networks and the internet.
- Overlay multiple graphs for comparisons
- Tabulate up to 99 tests results for analysis
- Previous tests are always available for review.
- Test platforms can be operated remotely.
- The ChemInstruments EZ lab software operates with ChemInstruments AR-1000,TT-1100,BP-1000,PT-1000, LT-1000,HSU-1000 and is included with HSR-1000 & PMA-10000.



System Requirements

- Personal Computer
- 64 MB of RAM
- Pentium processor, 200 MHZ or higher

- RS-232 comm. Port
- Hard Drive space...3MB